

Master Test Plan: ESP32-Based ECG AI System

Executive Summary

This report presents a comprehensive testing framework for an ESP32-based medical ECG monitoring device with embedded AI. We detail the required hardware, firmware, and AI model specifications; design exhaustive test scenarios covering physiological and device-failure conditions; define performance, accuracy, robustness, and usability metrics; outline regulatory and ethical requirements; and propose a scoring rubric. The goal is to ensure **safety, reliability and clinical effectiveness**. Key requirements include at least 500 Hz ECG sampling, $\leq 30 \mu\text{V}$ noise and $\geq 90 \text{ dB}$ common-mode rejection ¹, and a validated AI model (e.g. CNN) that meets sensitivity/specificity targets. The plan prescribes strict pass/fail thresholds (e.g. sensitivity $\geq 90\%$, latency $< 100 \text{ ms}$) with remediation steps. Reporting will include summary tables, workflow diagrams, and charts of expected metrics.

System Specifications

- **Hardware:** The device uses an **ESP32-WROOM-32** module (Dual-core 240 MHz, up to 8 MB Flash, up to 8 MB PSRAM) ². It must include an ECG analog front-end (e.g. **AD8232**) with rails-to-rail output and built-in filtering ³. The AD8232 (or similar) requires a 2.0–3.5V supply, draws $\sim 170 \mu\text{A}$, and provides $100\times$ gain with 80 dB CMRR ³, outputting a conditioned analog ECG signal. This output connects to an ESP32 ADC pin (e.g. VP) ⁴. Electrodes (2 or 3, adhesive chest patches) interface via leads with lead-off detection. Power is provided by a battery (spec TBD), with voltage regulation to 3.3 V. **Connectivity** (Wi-Fi/Bluetooth) must be tested for reliability and secure data transfer. Expected hardware constraints (e.g. ADC 12-bit resolution, $\pm 6\%$ error ⁵) should be documented as variables in the test prompt.
- **Firmware:** Real-time firmware must run on ESP-IDF or Arduino. Core functions include: continuous high-rate ADC sampling (target $\geq 500\text{--}1000 \text{ Hz}$ ¹), digital filtering (e.g. notch, high-pass for baseline removal), R-R interval computation, and data buffering. It must support wireless data uplink (with fallback if offline), power management (sleep modes), and user interface (LEDs, buzzer, or display). The firmware should implement safety checks (lead-off detection, low-battery alerts) and logging of events. Specify any RTOS or task scheduling (e.g. FreeRTOS tasks for sampling vs. networking). Firmware versioning and over-the-air update mechanisms should follow IEC 62304 processes ⁶.
- **AI Model:** The AI component (e.g. arrhythmia classifier) may run onboard (tinyML) or on connected app/cloud. Unknowns (model architecture, input window length, classes) are **parameters** in the master prompt. For example, *CNN_ArrhythmiaNet* might take X seconds of single-lead ECG and output rhythm class (AFib, PVC, etc). Model requirements: trained on a diverse ECG dataset (single-lead), validated on hold-out set. It should achieve high sensitivity/specificity for targeted rhythms (in Rajpurkar et al., a 34-layer CNN exceeded cardiologist sensitivity/PPV on single-lead ECG ⁷).

Assume metrics like ROC-AUC >0.95 for key classes as goals. Specify model input/output formats, expected latency, and memory footprint (it must fit on ESP32 if on-device inference).

- **Software/UW:** The system may include a user interface (mobile app or device display). UI specs: display real-time ECG waveform, heart rate, alerts. Controls for threshold settings, patient ID, start/stop. Logging of events. Interface must comply with usability guidelines (IEC 62366) ⁸ and privacy (HIPAA/GDPR).

Comprehensive Test Cases

Physiological Conditions

Define test ECG waveforms covering *normal and pathological rhythms*. For each, simulate at least 5–10 minutes of data and evaluate system output. Use an ECG simulator (hardware or software, e.g. PhysioNet ECGSYN or a Fluke ProSim) to generate waveforms. Examples:

- **Normal Sinus Rhythm:** at rest (60–100 bpm), deep breathing (baseline wander), lying vs. standing (heart rate changes).
- **Sinus Bradycardia/Tachycardia:** (~40 bpm, ~150 bpm).
- **Premature Beats:** frequent PVCs/PACs (single and paired), couplets, bigeminy, trigeminy.
- **Atrial Arrhythmias:** Atrial fibrillation (irregularly irregular, no P-waves), atrial flutter (sawtooth waves, 2:1/3:1 conduction), atrial tachycardia.
- **Supraventricular Tachycardia:** AVNRT/AVRT patterns (narrow QRS, sudden onset/offset).
- **AV Blocks:** 1st, 2nd (Mobitz I/II), 3rd degree (complete block with ventricular escape).
- **Ventricular Arrhythmias:** Ventricular tachycardia (monomorphic, polymorphic/Torsades), Ventricular fibrillation (coarse, fine VFib simulation for alarm testing).
- **Pacemaker/Artifact:** Paced rhythms (spikes on ECG), check if algorithm misclassifies. Also test benign abnormalities (ST-elevation simulating MI vs. electrolyte changes).
- **Hypotheticals:** Electrolyte imbalance effect (if supported), hyper/hypothermia wave changes.
- **Signal Variations:** Vary electrode placements (simulate chest vs. limb leads if applicable).

Motion and Artifacts

Test robustness to motion and noise:

- **Muscle (EMG) noise:** Add EMG-like interference (~20–200 Hz) at varying amplitudes.
- **Baseline Wander:** Simulate breathing artifacts, patient movement.
- **50/60 Hz Interference:** Intentionally expose to mains hum; verify CMRR performance (should suppress it ~80–90 dB ³ ⁹).
- **Lead Vibrations:** Mechanically jiggle electrode or connector to induce spikes/dropouts.

Electrode and Lead Faults

- **Leads-Off:** Simulate a disconnected electrode. AD8232 has leads-off detection; ensure firmware alerts and stops processing.
- **High Impedance:** Increase skin-electrode impedance (resistor network) to simulate poor contact; check signal degradation.

- **Short Circuits:** Short two leads (simulate cross-talk); device should warn or fail gracefully.

Power and Connectivity Failures

- **Battery Low:** Drop supply below threshold; verify battery monitor and system behavior (graceful shutdown or low-power mode).
- **Power Loss/Reset:** Remove power abruptly; test data integrity on reboot (log persistence).
- **Wireless Disconnect:** If using Wi-Fi/BLE, simulate range loss or reboot of access point; check auto-reconnect or local caching of data.
- **Firmware Crash:** Force a code exception; verify watchdog reset and data recovery.

Edge and Rare Cases

- **Extreme Physiology:** Use ECG of athletes (very high stroke volume) or infants (fast rates).
- **Comorbid Signals:** Feed ECG from a patient with other devices (like ICD spikes).
- **Data Anomalies:** Truncate data, insert NaNs; system should handle (reject, error-code).
- **AI Limits:** Provide ECG with unknown rhythms; check if model outputs "unclassified" or default.

These tests should be **prioritized by risk**: e.g., life-threatening arrhythmias (VTach/VFib) and sensor disconnections are top priority, followed by routine rhythms and user errors.

Performance and Accuracy Metrics

- **Sampling & Latency:** Verify ADC sampling ≥ 500 sps (target 500–1000 Hz) ¹. Measure end-to-end latency: time from ECG event (e.g. R-wave) to AI classification output. Target <100 ms for real-time alerts. If Wi-Fi enabled, factor in network delay.
- **Noise and Resolution:** Measure system noise floor (RMS or peak-to-peak). It should be $\ll 30 \mu\text{V}$ ¹. Confirm ADC resolution: ESP32 12-bit ADC yields $\sim 3.3 \text{ V}/4096 \approx 800 \mu\text{V LSB}$ ¹⁰ ¹¹ (limit test precision).
- **Sensitivity/Specificity:** For each arrhythmia class, compute sensitivity, specificity, PPV, NPV, and ROC-AUC. E.g. for AF detection, require sensitivity $\geq 90\%$, specificity $\geq 90\%$. (Rajpurkar et al. achieved >0.97 AUC on several classes ⁷.) Use statistically significant sample size (e.g. >100 episodes per class) to compute 95% confidence intervals.
- **Accuracy under Noise:** Add Gaussian noise and motion artifacts; measure drop in sensitivity/precision. Define *noise tolerance* (e.g. maintain $\geq 80\%$ sensitivity with X dB added noise).
- **Signal Delay and Skew:** If multiple leads (unlikely on single-lead device), ensure inter-channel skew $<100 \mu\text{s}$ ¹. For single-lead, ensure consistent timestamping of events.
- **Throughput:** If streaming data, measure maximum sustainable data rate (kbps) and CPU load. Ensure no buffer overflow at peak workloads.
- **Usability:** Rate UI on clarity and response time. Metrics: time-to-configure (should be <2 min), error rate (should be none in guided test), user satisfaction (qualitative).
- **Regulatory Compliance:** Check checklist for IEC 60601-1 (safety, leakage, EMC), IEC 60601-2-47 (ambulatory ECG standards), IEC 62304 (SW lifecycle), ISO 13485 (QMS), ISO 14971 (risk) ¹² ¹³ ¹⁴. For FDA/CE, verify needed submissions (e.g. 510(k) for ECG machine ¹⁵, MDR Class IIa/IIb per MDR Rule 11 ¹⁶ ¹⁷).

All measurement methods must follow IEEE/IEC test procedures where possible (e.g. IEC 60601-2-25 signal tests ¹). Calibration equipment (calibrated signal generator, multimeter) must be traceable.

Regulatory, Safety, and Accountability

- **Safety & Standards:** The device must meet IEC 60601 general and ECG-specific standards (60601-1, -2-25, -2-27, -2-47 ¹²). Examples: input impedance $\geq 2.5\text{ M}\Omega$, passband 0.05–150 Hz, CMRR $\geq 90\text{ dB}$, channel crosstalk $<2\%$ ¹⁰ . Documentation must reference these. Ensure compliance with ISO 13485 QMS ¹³ and ISO 14971 risk management throughout design ¹⁴ .
- **Software Lifecycle:** Follow IEC 62304 processes for software development and testing ¹⁸ ⁶ . Classify software safety level (A/B/C) based on risk (IEC 62304 safety classes). Maintain design history, version control, and verification evidence (unit tests, integration tests, validation).
- **Usability:** Perform human factors analysis per IEC 62366-1 ⁸ . Test the UI with representative users to catch use-error risks (e.g. misinterpreting alerts).
- **Cybersecurity & Privacy:** Ensure data encryption in transit and at rest. Follow IEC 81001-5-1 (cybersecurity lifecycle) guidelines ¹⁹ . Implement user authentication, secure pairing. Confirm GDPR/HIPAA compliance for personal health data (anonymize test records).
- **Traceability and Logging:** The system must log all events and data (raw ECG, classifications, errors) with timestamps. This audit trail aids post-market surveillance and investigation of incidents. Maintain test logs for FDA/Notified Body review.

Failure Modes & Mitigation: List potential failures (e.g. false negatives in arrhythmia detection, hardware battery failure, lost connectivity). For each, define mitigation: e.g. if heart rhythm not classified, default to safe alarm and advise manual ECG; if battery low, alert user early; if sensor off, halt alarms and log. Use ISO 14971 analysis to manage residual risks ¹⁴ .

Scoring Rubric (0–10)

Criterion	0–3 (Poor)	4–6 (Moderate)	7–10 (Excellent)
Accuracy	Many missed/false alarms; sensitivity $<80\%$.	Some detection, sensitivity 80–90%.	Near-perfect detection; sensitivity $\geq 90\%$, high PPV (e.g. $>90\%$). (Benchmarks per ⁷ .)
Latency	$>300\text{ ms}$ end-to-end delay.	100–300 ms (accept as borderline for real-time).	$<100\text{ ms}$ (real-time, negligible lag).
Usability (UI)	Confusing UI, high error rate, $>5\text{ min}$ setup.	Basic UI, some guidance needed, 2–5 min setup.	Intuitive interface, $<2\text{ min}$ setup, clear alerts.
Ease of Use	Poor documentation/training; users struggle.	Adequate docs; moderate training needed.	Clear docs/tutorial; self-explanatory use.
Innovation	Conventional design; no novel features.	Some innovation (e.g. unique sensor or algorithm).	Novel AI approach or breakthrough integration.

Criterion	0–3 (Poor)	4–6 (Moderate)	7–10 (Excellent)
Real-World Feasibility	Not tested in target environments; heavy hardware.	Laboratory-only demo; incomplete robustness.	Portable, power-efficient, validated in field.
Safety/ Accountability	No risk analysis; poor logging.	Basic risk checklist; partial logs.	Full ISO14971-based risk mgmt; complete audit trail.
Regulatory Preparedness	No plan for FDA/CE.	Aware of regulations but no evidence of compliance.	Formal strategy for 510(k)/MDR, demo compliance with IEC standards.

(Scores are averaged for overall “readiness”. Each criterion’s thresholds should be adjusted to project needs.)

Master Testing Prompt

Agent Prompt:

“You are an expert medical device test agent. First, list all required hardware and firmware specifications for the ESP32 ECG AI system (e.g. ADC bits, sampling freq, sensor type, connectivity, power). Then enumerate needed AI model details (architecture, inputs, outputs, metrics to measure). Next, design comprehensive test cases covering: physiological ECG conditions (normal rhythm, arrhythmias: AFib, VTach, PVC, blocks, etc.), motion artifacts (EMG noise, movement), electrode faults (lead-off, short), and system failures (battery low, connectivity loss). For each case, define expected outputs and metrics. Specify real-time metrics (latency, throughput), accuracy metrics (sensitivity, specificity, PPV, NPV, AUC), and robustness metrics (noise tolerance, drift). Include usability and UI tests (interface clarity, response time) and regulatory checks (IEC 60601-2-25/-47 compliance, ISO 13485/14971, FDA 510(k)/EU MDR classification ¹² ¹³). Address ethical/safety failure modes (e.g. false negative alarm) and mitigation. Provide a structured output: prioritized test checklist, equipment/phantom profiles, data analysis plan, pass/fail criteria with remediation, and reporting templates.”

This master prompt guides the agent to cover every aspect of evaluation (specs, tests, metrics, standards) in a systematic way, ensuring no critical area is omitted.

Prioritized Test Checklist and Measurements

1. Signal Quality & Safety:

- Verify ECG front-end gain, bandwidth (0.05–150 Hz) and noise $\leq 30 \mu\text{V}$ ¹ .
- Check electrode input impedance ($\geq 2.5 \text{ M}\Omega$) ⁹ , leads-off detection.
- Test CMRR with 50/60 Hz interference; expect $\geq 80 \text{ dB}$.

- Ensure patient leakage currents meet IEC 60601-1 safety limits.

6. Physiological Accuracy:

- Normal sinus waves: R-wave timing accuracy, heart rate error $< \pm 2 \text{ bpm}$.
- Arrhythmia detection: sensitivity/specificity for AF, VT, PVC, etc. (target $\geq 90\%$).

- 9. Rate detection: Tachy/brady detection thresholds and alerts.
- 10. Threshold crossing: correct identification of blocks, flutter.
- 11. **Artifact Robustness:**
 - 12. Baseline wander (amplitude % distortion), muscle noise (SNR degradation).
 - 13. Motion (walking, jogging): false-positive alarm rate $\leq X\%$.
 - 14. Electrode detach: system pause/sound alarm within 2 sec.
- 15. **Performance:**
 - 16. Latency: R-peak to classification <100 ms.
 - 17. Sampling: stable at 500–1000 Hz during Wi-Fi/Bluetooth use ²⁰.
 - 18. Data throughput: sustained without loss for continuous monitoring.
- 19. **Power/Connectivity:**
 - 20. Battery endurance: at least Y hours of continuous use.
 - 21. Low battery alert at 20% remaining.
 - 22. Reconnect after wireless dropout within 30 s.
- 23. **Usability:**
 - 24. UI: legibility (font size ≥ 12 pt), color contrast, touch response ≤ 200 ms.
 - 25. Menu navigation: complete key tasks in <2 steps.
 - 26. Instruction clarity: 100% of test users can start monitoring after 5 min tutorial.
- 27. **Regulatory & Documentation:**
 - 28. Risk analysis: Hazards/risk mitigation documented (ISO 14971).
 - 29. Compliance test: Reference IEC/FDA test methods (e.g. IEC 60601-2-25 Figure 201.110 waveform reproduction ²¹).
 - 30. Audit trail: Logging enabled (event, time, ECG data).
- 31. **Failure Modes:**
 - 32. Simulate wrong classification: verify safe fallback (e.g. persistent alarm requiring review).
 - 33. Data corruption: test integrity checks, no undetected crash.

Required Equipment and Patient Profiles

- **ECG Signal Generator:** Calibrated waveform simulator (e.g. Fluke ProSim, Biopac) capable of generating AAHA/ANSI-standard ECG rhythms and noise. If unavailable, use software (PhysioNet ECGSYN ²²) to produce digital signals fed through a high-speed DAC to the ADC input.
- **Electronic Loads:** To simulate electrode impedance changes and shorts (resistor/capacitor network per IEC).
- **Phantom/Simulator:** If available, use a patient simulator (NETECH, Philips) for holistic vitals (ECG, BP, etc.).
- **Data Acquisition:** High-resolution scope or DAQ (≥ 16 -bit, 1000 Hz) to capture raw ECG for ground truth.
- **Computing:** PC with analysis software (MATLAB/Python) to run analyses (signal metrics, ROC curves).

Simulated Patient Profiles:

- **Healthy Adult:** Sinus rhythm at 60 bpm, moderate noise (simulate relaxation vs. exercise).
- **Elderly Arrhythmia Patient:** Frequent AFib and PVCs; slower conduction.
- **Athlete:** Bradycardia (40 bpm), high QRS amplitude.
- **Pediatric:** High rate (120 bpm), small amplitude.
- **Pacemaker Patient:** Paced spikes present (use simulator mode).
- **Artifact-Prone:** Simulate profuse muscle noise (exercise), sweaty skin (variable impedance).

Each profile should run multiple sessions (rest vs. activity) to collect robust data.

Data Collection and Analysis Plan

- **Collection:** For each test, log raw ECG input and device outputs (heart rate, arrhythmia label, alarms). Synchronize timestamps. Store logs securely.
- **Repetition:** Repeat key tests ($n \geq 3$) to account for variability.
- **Analysis:** Use statistical scripts to compute mean latency, error bars for metrics. Plot confusion matrices for classification. Generate ROC curves and compute AUC for each rhythm.
- **Visualization:** Create tables of sensitivity/specificity per test, time-to-event charts (delay from arrhythmia onset to alarm). Plot distribution histograms of noise, latency (target bell-shaped around ideal). (Example chart: superimposed ROC curves for AFib vs. VTach.)

All code and data must be version-controlled. Ensure traceability from raw data to final metrics.

Pass/Fail Thresholds and Remediation

- **Noise/Resolution:** PASS if system noise $\leq 30 \mu\text{Vpp}$ and ADC LSB $\leq 5 \mu\text{V}$ equivalent (ideally); otherwise: improve filtering or ADC oversampling.
- **Sampling Rate:** PASS if ≥ 500 Hz stable; else: reduce Wi-Fi load or lower ADC attenuation (per ²⁰).
- **Accuracy:** PASS for each arrhythmia if sensitivity $\geq 90\%$ and specificity $\geq 90\%$. If FAIL, consider retraining model, adding more data, or adjusting decision thresholds.
- **Latency:** PASS if < 100 ms. If 100–300 ms, optimize code (e.g. use DMA, simplify model) or relax requirements if clinically acceptable. > 300 ms is FAIL (increase CPU or offload AI).


- **Connectivity:** PASS if device auto-reconnects within 30 s and caches <1 min of data. FAIL: implement buffering and retry logic.
- **Battery:** PASS if operational >8 hours on design battery. If <6 h, investigate power profile, enable deeper sleep modes, or increase battery capacity.
- **Usability:** PASS if $\geq 80\%$ of test users complete key tasks without error. If usability issues arise, revise UI (follow IEC 62366 methods ²³).

Remediation steps should be documented in a corrective action report, feeding back into development.

Reporting Template

The final report should follow this structure:

- **Executive Summary:** Key findings, overall scores (0–10) for each criterion (see rubric), major risks identified, and go/no-go recommendation.
- **System Specifications:** Itemized hardware/firmware/AI model specs used in tests, with citation of standards and datasheets (e.g. ESP32 ADC, AD8232 features ³ ²⁰).
- **Methods:** Description of test setup, equipment, signal sources, and participant profiles (simulate demographics and arrhythmias). Include a **Mermaid** workflow diagram of the testing process:



- **Results:** Use tables and charts. For example, a table comparing **performance metrics** across test cases:

Condition	Sensitivity	Specificity	Latency (ms)	Outcome
Normal Sinus	98%	99%	50	PASS
AFib (onset)	92%	95%	80	PASS
Ventricular Tachycardia	90%	93%	120	PASS (border)
Lead Off	N/A	N/A	Detected in 2s	PASS
Battery Low	N/A	N/A	Alert after 5h	PASS
etc.				

Include a **Mermaid timeline** showing test execution over calendar days:


```

gantt
  title Test Plan Timeline
  dateFormat YYYY-MM-DD
  section Planning
    Define Requirements      :done, des1, 2026-04-01, 3d
    Procure Equipment       :done, des2, 2026-04-04, 7d
  section Execution
    Physiological Tests     :active, ext1, 2026-04-12, 10d
    Failure Mode Tests      :ext2, after ext1, 7d
  section Analysis
    Data Analysis           :2026-04-22, 5d
  section Reporting
    Draft Report            :2026-04-27, 4d
    Review & Finalize       :2026-05-01, 2d

```

And a **chart** (e.g. a ROC curve or histogram) illustrating expected metric distribution (e.g. distribution of latency under load).

- **Findings:** Textual analysis of each test area, citing any deviations. For example, “ADC noise was measured at 50 μ Vpp (above spec), likely due to ESP32 internal noise ¹ – recommend adding external low-pass filter.”
- **Pass/Fail Matrix:** Summarize each criterion with PASS/FAIL and rationale.
- **Scoring:** Fill in scores (0–10) per rubric section with brief justification (e.g. “Accuracy 8: Sensitivity 92%/Specificity 94% (meets target). Latency 7: average 120 ms (target <100 ms).”).
- **Action Items:** List remediation needed for any FAILs or borderline results.

All sources (ESP32 specs, sensor datasheet, standards, and AI studies ³ ¹ ²⁴) must be cited when relevant. The executive summary should highlight readiness for clinical deployment or note critical issues.

¹ ⁹ ¹⁰ ¹² ²¹ A Practical Guide to ECG Device Performance Testing According to International Standards
<https://www.mdpi.com/2079-9292/14/19/3878>

² ESP32 Wi-Fi & Bluetooth Modules I Espressif
<https://www.espressif.com/en/products/modules/esp32>

³ AD8232 Datasheet and Product Info | Analog Devices
<https://www.analog.com/en/products/ad8232.html>

⁴ IoT-based Heart Monitoring with ESP32 | Espressif Systems
https://www.espressif.com/en/news/ESP32_Heart_Monitoring

⁵ ¹¹ ²⁰ Analog-to-Digital Converter (ADC) - - — ESP-FAQ latest documentation
<https://docs.espressif.com/projects/esp-faq/en/latest/software-framework/peripherals/adc.html>

⁶ ¹⁹ AI Device Standards You Must Know - ISO 13485, 14971, 62304
<https://www.hardianhealth.com/insights/regulatory-ai-medical-device-standards>

⁷ [1707.01836] Cardiologist-Level Arrhythmia Detection with Convolutional Neural Networks
<https://arxiv.org/abs/1707.01836>

8 23 IEC 62366-1:2015 - Medical devices — Part 1: Application of usability engineering to medical devices
<https://www.iso.org/standard/63179.html>

13 ISO - ISO 13485 — Medical devices
<https://www.iso.org/iso-13485-medical-devices.html>

14 ISO 14971:2019 - Medical devices — Application of risk management to medical devices
<https://www.iso.org/standard/72704.html>

15 21 CFR § 870.2340 - Electrocardiograph. | Electronic Code of Federal Regulations (e-CFR) | US Law | LII / Legal Information Institute
<https://www.law.cornell.edu/cfr/text/21/870.2340>

16 17 MDR Examples: Classification Examples for EU MDR Devices - OpenRegulatory
<https://openregulatory.com/articles/eu-mdr-classification-examples>

18 IEC 62304:2006 - Medical device software — Software life cycle processes
<https://www.iso.org/standard/38421.html>

22 ECGSYN - A realistic ECG waveform generator v1.0.0 - PhysioNet
<https://physionet.org/content/ecgsyn/>

24 Artificial Intelligence in Software as a Medical Device | FDA
<https://www.fda.gov/medical-devices/software-medical-device-samd/artificial-intelligence-software-medical-device>